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## **Microtracers – reliable checks on homogeneity and carry-over**

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COUNTING PARTS INSTEAD OF MEASURING CONCENTRATION

*The homogeneity and carry-over of animal feed can be checked with traditional tracers, such as manganese or cobalt, but it can also be measured with Microtracers. This latter method is based on particle counts and not on measuring the concentration.*

During the production of animal feed, it is extremely important that every batch has the same composition. At the same time, it is important to reduce the slippage of feed components to a minimum. For the animal-feed producer it is essential to check these criteria – if feeds are not mixed homogeneously, it costs money, and an incorrect composition of feed leads to inefficient animal production and to a waste of feed ingredients.

Traditionally, manganese or cobalt are used as tracers to measure carry-over, but since July 2006, Microtracers are also permitted (GMP+, 2006). The special aspect of microtracer technology is that it is based on particle counts and not on measuring a concentration. Cobalt, manganese, methyl violet and Microtracer RF are indirect external tracers, whereas Microtracer F, FS and FSS are direct external tracers.

### **Microtracers**

In Europe, Microtracer FSS and F and the services of Micro Tracers Services Europe (MTSE) are marketed by Jadis Additiva.

Microtracers consists of irregularly shaped iron particles. They can be supplied in a range of colours, which means that different feed components can be tracked in one mix line. The irregular shape of the particles allows the Microtracers to travel in the feed mix and prevents them from drop-

ping down during the production process.

Microtracers are applied in mixer tests to check uniformity and to optimise the mixing time and the addition time of the micro components. They can also establish the degree of carry-over and they can demonstrate feed additives – e.g. coccidiostatics, selenium, medicated feeds. Microtracers can also be used as a marker for measuring the passage speed through the intestines. Tracers need to be stable and must correspond with the smallest relevant particles in the mixture. The analysis of the tracer must be accurate, simple and economical, and ideally it should be possible to carry it out in situ.

Research by TNO in Zeist demonstrated that Microtracers produced reliable results and that they are simpler and more efficient in application than manganese and cobalt.

### **Mason Jar-method**

Microtracers are added to the mix line as a concentrate or by means of an automated microdosage device, and they are added at the point the other components are also added.

The amount of tracer must be equal to the amount of the least concentrated component of the mixture. For pre-mixes this is 1:100 000, for complete feeds it is 1:10 000. MTSE and Jadis Additiva have invested substantial sums in

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
the development and optimisation of their sampling and analysis methods.

Correct sampling is essential to an accurate analysis. The number of samples, the amount of sample, the time of sampling and the equipment used determine the analysis result to a considerable extent.

Two methods have been developed for the analysis of the number of particles – the ‘Mason Jar Technique’ and the ‘Rotation-detector technique’. MTSE provides support for the implementation of the detection technique in the form of workshops. The organisation also signs agreements with accredited laboratories in the Netherlands (STERLABS) and organises ring tests to safeguard the quality of the analyses.

In the Mason Jar method, the qualitative method, the feed sample (ground pellets) is weighed and transferred into a 500 ml jar with a magnetic lid. Filter paper is inserted into the lid and the jar is shaken for one minute. Then the filter paper is removed, and sprayed with an ethanol/water mixture. After 30 seconds the filter paper is removed from the lid and dried.

Now it is possible to count the coloured particles on the filter. This method is easily applied in the factory; it does not require special equipment.

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The Rotary detector filters a weighed sample; after which the iron particles can be demagnetised.

**Rotary Detector**

The Rotation-detector technology, the quantitative method, requires a rotation detector, a demagnetiser and filter paper (70mm and 240mm). The weighed sample is filtered in the rotation detector. The iron particles are demagnetised on the 70mm filter paper. The 240mm filter paper is sprayed with an ethanol/water mixture and the sample is transferred from the 70mm filter to the 240mm filter. The big filter is dried, and then the individual coloured particles are counted.

According to the EU Directive, the mixture accuracy must be 1:100 000 for pre-mixes. When adding 10 g of microtracer F (with 25 000 particles per gram) per tonne of feed, the feed contains 250 particles per kg. In a feed sample of 200 g there would be 50 particles. This quantity can be counted manually (table 1).

The statistical evaluation of the analysis results is substantiated properly by MTSE and Jadis Additiva.

As the method is based on a particle count (Poisson distribution), and not on a concentration measurement, the interpretation of the values requires different statistical methods.

**Übersicht 1: Mikrotracer**  
**Table 1: Microtracer**

Product/product	Teilchen pro Gramm/ particles per gram	Teilchengröße/ Particle Size
Mikrotracer F/ Microtracer F	25 000	150-300 mym
Mikrotracer FSS/ microtracer FSS	150 000	75-300 mym
Mikrotracer RF/ Microtracer RF	> 2 000 000	75-150 mym
Farben/Colours:	blau, rot, orange, grün, violett und Mischfarben	blue, red, orange, green, violet and mixing colors
Wasserlöslich – nicht wasserlöslich/ watersoluble – not watersoluble		
Hitzebeständig/heat resistant		
Analyseverfahren/ analysis method:	15 Teilchen/particles	qualitativ/qualitative
	16-100 Teilchen/ particles	quantitativ und manuell/ quantitative and manual
	101-500 Teilchen/ particles	quantitativ/ quantitative (Scanner/ scanner)
	> 500 Teilchen/ particles	quantitativ/quantitative (Waage/scale)

Quelle: Jadis